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*Lead the way, all the way!*

GENERAL WARNINGS

- Always maintain the trolley well, in particular checking tyre pressure, as this affects the performance of the brakes.
- When getting on or off the wheelchair, make sure that both brakes are engaged.
- If the trolley is not being pushed or used, make sure that both brakes are engaged.
- Care must be taken when negotiating kerbs and making uphill turns.
- Steep and uneven surfaces should be avoided as this can cause the trolley to tip over.
- Avoid soft surfaces as the wheels may collapse.
- Move the trolley smoothly in all conditions.
- Always observe your surroundings to identify hazards, such as furniture, doors, holes, hills, kerbs, entrances, etc.
- Keep both feet on the footrests at all times, and ensure that clothing is properly tucked and does not block the wheels.
- Always drive carefully on hills.
- Always keep both hands on the push handles for proper control of the trolley (applies to helper).
- Always remember to maintain balance.
- Check the frame and wheelchair components regularly and always follow the instructions given in this manual during and after the warranty period.
- Tie the upholstery or replace it with new upholstery if it becomes loose.
- Equip the trolley with anti-tipping devices to avoid falls backwards.
- The wheelchair should be fitted with safety belts designed for disabled people.
- The wheelchair should be fitted with heel loops for people with disabilities.
- Wash both hands after greasing the trolley.
- Keep the trolley away from sources of fire.
- Under no circumstances should the trolley be overloaded.

- Do not hang any heavy objects on the rear of the trolley, as this may affect the stability of the trolley especially when climbing hills.
- The wheelchair can only be used by one person.
- Do not use the wheelchair when the patient is unconscious.
- Do not use other means of transport to push the trolley.
- Do not lift the front wheels while driving, as there is a risk of falling and causing serious injury.
- Do not get into the trolley while it is moving.
- Do not move the trolley too fast and do not make sudden turns.
- Do not drive the trolley across a gradient (in either direction) greater than 7 degrees.
- Do not stand on the footrests when getting in or out of the wheelchair.
- Do not climb stairs using a wheelchair, either alone or with a carer.
- Do not attempt to lift the wheelchair with moving (detachable) parts. Lifting with the help of moving parts of the wheelchair may lead to injury to the user and/or carer or damage to the wheelchair.
- Do not use the pushchair if the backrest is folded. This will change the centre of gravity of the pushchair and may cause the pushchair to tip backwards. The backrest should only be folded when storing or transporting the pushchair.
- Do not hang any objects from the rear joint, as this may cause the backrest to fold down and cause serious injury to the user.

#### **KARMA AS A SERVICE PROVIDER**

- Tools are not required to assemble the trolley. The trolley should be assembled by an authorised Karma dealer. The authorised Karma dealer is also the service provider for the pram. Unauthorised repairs, upgrades and additions will void the warranty and may cause injury. Any adjustments to the rear wheel position should be carried out by an authorised Karma service provider. The service provider has the necessary tools for the wheelchair.
- If the wheelchair has damaged components or loose parts, contact your dealer for repair and/or replacement of the wheelchair. Always use parts authorized by KARMA. Using parts not authorized by KARMA may void the warranty.
- Customers should turn to their dealer for all maintenance services so that the warranty they have obtained is not voided. Authorised service providers should always be approached. Unauthorised repairs will invalidate the warranty. Procedures other than those described in this manual must be carried out by qualified technical personnel. These personnel are in possession of service manuals with a complete list of maintenance services required to keep the wheelchair in the best possible condition. Refer to the warranty section to find out what KARMA guarantees and what maintenance services and parts the user will have to pay for themselves.
- KARMA does not currently have an official list of service providers. Please contact your dealer. If you are unable to contact your dealer, please contact an authorised Karma dealer. If this proves impossible, please contact us and we will try to arrange a service repair. Please note that there will be a charge for collection and delivery of the trolley.

CE marking

This wheelchair complies with the Medical Device Directive 93/42/EEC.

**WARNING:**

Before purchasing and using a Karma trolley, please contact a qualified Karma technician to ensure correct trolley selection and to receive training on use. You may need to install additional equipment on the trolley or change the configuration. Please contact your Karma service provider.

**TARGET PATIENT GROUP**

Persons suffering from diseases, dysfunctions or injuries for the treatment/rehabilitation or compensation of which this device is intended (see section on intended use of the device in these instructions for use). The device can be purchased by the user independently or on the recommendation of a doctor, therapist or other specialist. Whether you purchase the device yourself or on the recommendation of a doctor/therapist/other specialist, you must take into account the available sizes/necessary functions and variants of the device, the indications and contraindications for use and the information provided by the manufacturer.

**1. INTRODUCTION**

Please read this user manual carefully before using the wheelchair. Incorrect use of the wheelchair may lead to health damage, injury or accidents. Thus, in order to use the wheelchair safely and successfully, please read this manual carefully.

- The user manual contains instructions on how to operate the wheelchair, assembly instructions, as well as instructions on how to prevent possible accidents.
- This user manual is intended for the wheelchair: Ergo Lite (KM-2501& KM-2512 ).
- The meaning of the individual graphic symbols used in the text is explained below. Pay particular attention to the parts of the text indicated by these symbols.

 WARNING:	Incorrect use can lead to serious injury or death.
 NOTES:	Incorrect use can lead to personal injury and/or damage to the wheelchair.
 SUGGESTION:	The instructions given must be followed in order to keep the trolley in good and proper condition.

- The user manual includes a repair and maintenance diagram, as well as a warranty.
- If there is another person using the wheelchair, make sure they are familiar with this user manual.
- As trolley designs change, some of the illustrations and photos shown in this manual may not correspond to the type of trolley purchased. We reserve the right to make modifications to the designs without prior notice.

**2. SAFETY**

**2.1 Entrances, doors**

The safest way to operate the trolley is on a level floor free of any obstacles. If possible, always use ramps when entering or leaving the building. Keeping your hands on the push handles allows maximum balance. Do not force yourself around obstacles by pushing and never use the sides of the door to push off.

**2.2 Leaning**

Do not lean over the side of the wheelchair as this can lead to loss of balance and falling. (See Fig. 2.3).

**2.3 Moving up hills**

- Do not make turns on hills because of the risk of overturning.
- When assisting, move straight up. If it is necessary to move downwards, walk backwards until the user is on level ground. (See Fig. 2.4 and 2.5).
- Avoid stopping on hills; this can lead to loss of control of the trolley.
- Do not use the brakes to reduce speed or stop; this can lead to locking of the brakes and loss of control of the wheelchair.

- If possible, move up hills together with a caregiver.
- Operating restrictions: do not climb hills of more than 5 degrees without the assistance of a caregiver.

## 2.4 Curbs

- If possible, a carer should always be used (if moving with a carer, the wheelchair must have push handles).
- Avoid strong shocks when stepping off the kerb. Shocks can have an adverse effect on the trolley.
- Under no circumstances should you attempt to negotiate a kerb whose height exceeds the normal height (30 - 50 cm).
- Always exercise caution when crossing kerbs and do not drive over a kerb in reverse without the assistance of a carer.

## 2.5 Stairs

- Always use two or more caregivers. Tilt the wheelchair to the balance point. One carer (at the back) holds the pushchair on the first step resting their arms firmly for lifting. The second carer, holding the fixed part of the front frame, lifts the wheelchair above the steps while the first carer places one foot on the next step and repeats the action. (See Fig. 2.1).
- Do not take the stairs going forward.
- Do not use escalators. Use a lift.

Fig. 1


Fig. 2

## 2.6 Streets

- Always take extra care when crossing the street. Reflective tape on the trolley and/or on clothing provides better visibility to drivers.
- Avoid holes and uneven surfaces that can cause the wheels of the wheelchair to lock. (See Fig. 2.6)

## 2.7 General information

- Always check that BOTH brakes are engaged when getting on or off the wheelchair.
- Do not stand on the footrests when getting in or out of the wheelchair.
- It is not permitted to carry more than one passenger. Karma wheelchairs are designed for one user.
- Please take care that the clothing does not block the wheels of the trolley.

## 2.8 Upholstery

- Use warm, soapy water to clean the upholstery.
- Upholstery cleaner/foam can also be used, but avoid using solvents.

## 2.9 Anti-tipping devices

Please purchase anti-tipping devices and always use them to avoid injury. The distance between the anti-tip device and the wheels is more than 5 cm; the floor height is not more than 5 cm (See Fig. 2.2).

## 2.10 Maximum user weight limit

Please refer to the section "Maximum User Weight" in Section 8 SPECIFICATIONS. Loads in excess of the maximum load capacity of the trolley may damage the trolley and cause faults leading to serious danger. The warranty does not cover damage caused by incorrect operation of the trolley.

## 2.11 Avoid misuse

- Please be aware that any vehicle can cause injury if operated incorrectly. This manual contains a number of tips to enhance user safety. Incorrect use of the wheelchair is dangerous to both the user and bystanders. ALWAYS obey traffic regulations and the instructions given in this manual. The wheelchair can serve the user well for a long time if used correctly.
- Do not allow children to stand or play on the trolley. (See Fig. 2.7).

Fig. 1

Fig. 2

Fig. 3


Fig. 4

Fig. 5

## 2.12 Labelling

Before using the trolley, please read the labels on the trolley carefully. Please do not remove the labels and keep them for future reference.


3. PURPOSE AND COMPONENTS

3.1 Purpose


- Ergo Lite (KM-2501 & KM-2512) is a transit wheelchair. This means that it can only be propelled by the carer (transit).
- The wheelchairs are designed for users weighing up to 100 kg. The rear wheels are equipped with brakes for the carer mounted on both sides.
- The pushchair is equipped with pneumatic tyres for good traction and comfort. The pushchair also has an extremely sturdy aluminium frame with a set of accessories and cannot withstand extremely high temperatures (neither the user nor the carer need to touch metal parts when using the pushchair).
- The front wheels are fitted with PU tyres for extra durability; pneumatic front wheels are available on request.
- When prescribing a wheelchair for use by a person with a partial or full amputation (above or below the knee, one or two limbs) or other conditions that compromise the natural centre of gravity, we recommend that a therapist or dealer is contacted and that anti-tip devices are fitted.

### 3.2 Description of the KM- 2501 trolley parts

- 1 Push handle
- 2 Keeper brake
- 3 Rear connector
- 4 Rear wheel
- 5 Anti-tipping device (Option)
- 6 Transverse locking bar
- 7 Front wheel
- 8 Backrest upholstery
- 9 Armrest
- 10 Side panel
- 11 Seat upholstery
- 12 Calf support
- 13 Footstool

### 3.2 Description of the KM- 2512 trolley parts

1. push handle
- 2 Keeper brake
3. rear connector
4. rear wheel

5. Anti-tipping device (Option)
6. transverse locking beam
7. front wheel
- 8 Backrest upholstery
9. armrest
10. side panel
11. seat upholstery
12. Footrest
13. calf support
14. foot rest

#### 4. OPERATION

##### 4.1 Preparation of the trolley

##### 4.1.1 Disassembly

- Holding the armrests, pull them away from you (See Fig. 4.1).
- Pointing your fingers inwards while keeping your hand on the frame, push downwards until the seat lowers completely. Make sure the frame is securely seated in the handle. (See Fig. 4.2).
- Turn over the footrests (see Fig. 4.3).
- The trolley is ready for use. (See Fig. 4.4).

Alternatively, you can stand at the side of the wheelchair and carry out step 1 by pulling the nearest armrest towards you and pushing the other armrest in the opposite direction.

Fig. 4.1

Fig. 4.2

Fig. 4.3

Fig. 4.4

#### **WARNING:**

When performing step 2, do not put your fingers around the sides or under the seat as they may be pinched between the seat and the frame. (See Fig. 4.5)


Fig. 4.5

**WARNING:**

The seat upholstery must be fitted before using the wheelchair. (See Fig. 4.6)

Fig. 4.6

#### 4.1.2 Folding

- Fold down the footrests (see Fig. 4.7)
- Remove the upholstery by holding it in the middle in the front and rear sections (See Fig. 4.8)
- Press the two release levers on the backrest frame to fold the backrest for storage or easy transport. (See Fig. 4.9)
- Folding is now complete. (See Fig. 4.10) (The trolley cover is an optional accessory).

Fig. 4.7  
4.10

Fig. 4.8

Fig. 4.9

Fig.

**WARNING:**

Take care not to press your fingers together when folding the trolley.

#### 4.1.3 Regulations

##### 4.1.3.1 Rear wheel

The trolley is equipped with wheels for quick disassembly. If necessary, follow the instructions below.

##### 4.1.3.1.1 Removing the rear wheel

- Hold the trolley by the push handle.
- Press the hub cap and hold the wheel release mechanism.
- Remove the wheel. (See Fig. 4.11)

#### 4.1.3.1.2 Rear wheel assembly

- Hold the wheel. Holding the centre of the rear wheel by two adjacent spokes using your thumb, press the release mechanism. Take the axle of the quick release mechanism and press the button in the centre. This will cause the bearing to move to the opposite end of the axle. This will allow the axle to move through the centre of the wheel hub. The wheels are now ready to be fitted to the trolley frame.
- Check that the trolley brake is in the OFF position. Take the assembled wheel and axle and press the button in the centre of the axle as above. This will cause the bearings to come apart and allow the axle to be inserted into the wheel clamp, which is located at the bottom of the wheelchair's rear frame. When you remove your thumb from the button, the bearings should protrude from the axle at the rear of the wheel clamp. This further secures the wheel in the frame. (See Fig. 4.12)

#### WARNING:

Ensure that the axle is fully in the wheel clamp of the trolley wheel and that the ball bearings are fully open and correctly fixing the wheel in the trolley frame.


Fig. 4.11


Fig. 4.12

#### 4.1.3.1.3 Recommended and unrecommended activities

- Keep axles clean and lubricated with oil or maintenance fluid.
- Ensure correct tyre pressure if pneumatic tyres are used.
- Do not attempt to remove the wheels if there is a user on the wheelchair and the wheelchair is fitted with anti-tip devices.

#### 4.1.3.2 Footrest height

- To adjust the height of the footrest, use the Karma tool to loosen the screw on the front of the footrest arm (See Fig. 4.13).
- Adjust the footrest to the required height by tightening the screw (see Fig. 4.13)

Fig. 4.13

Fig. 4.14

#### SUGGESTION:

- Ensure that the height of the footrest is adjusted correctly according to the length of the user's calf to reduce pressure on the body.
- The height of the footrest must be greater than 5 cm (see Fig. 4.14)

#### 4.1.4 Transport and storage

##### **SUGGESTION:**

Please store the trolley in a place that is not exposed to direct sunlight, rain and dew. For more information, please consult your Karma dealer.

##### **NOTES:**

Do not place any objects on the trolley when storing the trolley. Placing objects on the folded or unfolded trolley will cause damage that is not covered by the warranty. (See illustration below)

#### 4.2 User transfer

- Apply both handbrakes.
- Turn the footrests upwards and twist to the side.
- The user lowers onto the seat using the armrests as a backrest.
- Turn the footrests downwards so that the user can place their feet on them. Check that the heels are resting on the heel rests and not touching the front wheels.
- Release both handbrakes.

##### **WARNING:**

Always check that the front wheels are aligned with the rear wheels and pointing forwards, and that the rear of the wheelchair is leaning against a wall or held securely by the carer.

#### 4.2.2 Transfer from trolley

- Ensure that the trolley is close to the destination point.
- Withdraw the wheelchair appropriately so that the front wheels are pointing forwards, i.e. in line with the rear wheels. Ensure that the rear of the wheelchair is leaning against a wall or, if this is not possible, is held by the carer.
- Apply both handbrakes.
- Turn the footrests upwards.
- Move the user to the edge of the seat.
- Place one foot on the ground with the knee at a right angle and the other foot slightly to the side and behind.
- Place your hands on the front of the armrests and, shifting your body weight forward, move your head and shoulders above your knees pushing upwards.

##### **WARNING:**

- Do not stand on the footrests when getting in or out of the wheelchair (See Fig. 4.15).

- Use the carer's brakes when getting in or out of the wheelchair.
- Do not allow the tread on the tyres or the brake pads to become worn. If pneumatic tyres are used, they must be inflated to the correct pressure.

When transferring the user to or from the wheelchair, ensure that the brakes are engaged so that the wheelchair does not move. Also check that the footrests are turned upwards to make it easier to get on or off. (See Fig. 4.16)

### 4.3 General use

#### 4.3.1 Pushing the wheelchair by the carer

- Always check that the brakes are applied when the wheelchair is not being pushed by the carer. (See Fig. 4.17)
- Always check that the user's feet are freely placed on the footrests and that clothing is not blocking the wheels of the wheelchair. (See Fig. 4.18)
- When pushing, always use both hands to have good control of the trolley. (See Fig. 4.19)
- Steep or uneven surfaces should always be avoided as the trolley will tip over (See Fig. 4.20)
- Avoid excessively soft surfaces as the front wheels of the wheelchair may lock (See Fig. 4.21)
- Do not move the wheelchair too fast or make sudden turns - try to move smoothly.
- Stop completely and look both ways before crossing the street.
- Railway tracks should be crossed perpendicularly so that the wheels do not lock into the rails.
- Do not run the trolley at full speed on railway rails.

PLEASE have full contact with the person in the wheelchair. If stopping in front of a display window or when talking to another person, please remember to turn the wheelchair so that the wheelchair user does not have to turn on their own.

#### NOTES:

Do not lock the brakes too tightly as this can cause the user to tip over.

Fig. 4.15

Fig. 4.16

Fig. 4.17

Fig. 4.18

Fig. 4.19

Fig. 4.20

Fig. 4.21

Fig. 4.22


Fig. 4.23

#### 4.3.2 Brake

Please note: The 3-Step-Lock brake is a multifunctional brake that has the function of reducing speed and parking the wheelchair. Please read the brake instructions before using the brake to ensure your comfort and safety.

##### 4.3.2.1 Decreasing trolley speed

When moving up a hill, the carer may need to reduce the speed of the wheelchair to ensure safety. The first and second stage of the brake lock allow the speed of the wheelchair to be reduced.

To activate the first brake lock step, press the lever (see Fig. 4.23) until a single click is heard. To activate the second brake lock stage, press the lever (see Fig. 4.23) until two clicks can be heard.

##### 4.3.2.2 Parking the trolley

The third brake lock stage provides the parking function of the wheelchair. If the brake is not in any locking position, to activate the third brake locking stage press the brake lever until three clicks can be heard.

If the brake is in the first or second locking position, press the levers until you hear two or one click.

**Note:** If the tyres have worn out after a period of use, the brake must be adjusted or the parts replaced to ensure that the wheelchair parks correctly. Please contact your Karma dealer for assistance.

##### 4.3.2.3 Brake release

Regardless of the degree of brake lock, the brake release lever must be pressed to release the brake (see Fig. 4.23)

#### Attention:

3. Always check the brake to ensure it is fit for purpose.
4. Check the function of the brake release lever and the brake itself from time to time. When the lever or brake starts to fail, contact your Karma dealer for repair.

#### 4.3.3 Tilting the wheelchair by the carer

- Ensure that both of the user's feet are on the footrests and do not slip off. Ensure that the hand grips are secure and do not slip.
- Release both brakes.
- Grip the push handles firmly.
- Place a firmer foot on the appropriate tilt lever at the rear of the wheelchair.

- Notify the person in the wheelchair to tilt the wheelchair backwards.
- Pull the pushchair backwards with your hands and push downwards with your foot. This will allow the trolley to rest on the rear wheels. (See Fig. 4.22)

**WARNING:**

- When the wheelchair is inclined by about 25 degrees it will keep its balance on its rear wheels and be easily controlled. If it is inclined even more it will be difficult to keep its balance and the person in the wheelchair will not feel safe.
- Lower onto an upright back; keep your foot on the anti-tip device and slowly lower the trolley to the ground.

#### **4.3.4 Stepping off the kerb**

Note: It is best to avoid kerbs by using kerbs with ramps. This activity is best done with the help of a carer.

##### **4.3.4.1 Sliding forward from kerbs with the help of a carer**

**WARNING:**

- This is a more complicated method and may shock the user, but it allows better observation of traffic.
- Take care that the anti-tip devices do not get blocked on the kerb causing a fall.

#### **DO NOT ALLOW THE WHEELCHAIR TO TILT FORWARD WHEN STEPPING OFF THE KERB.**

- Position the wheelchair so that the front wheels are close to the edge of the kerb and ensure that there is no traffic on the road.
- Move the wheelchair slowly forward while tilting it slightly backwards. Be prepared to take the weight of the pushchair when the rear wheels of the pushchair cross the edge of the kerb. Ensure that both rear wheels of the pushchair touch the ground at the same time to avoid the pushchair falling.
- Gently lower the front wheels. **AVOID SHAKING THE TROLLEY.**

##### **4.3.4.2 Backwards stepping off the kerb with assistance from the carer**

**WARNING:**

This is an easier and more convenient method, BUT requires the carer to step backwards onto the carriageway.

- Back the trolley out to the kerb and ensure that both wheels are an equal distance from the kerb.
- Holding the push handles, slowly lower the rear wheels of the pushchair onto the roadway. Be prepared to control the weight of the pushchair when lowering and ensure that the rear wheels touch the ground at the same time.
- Tilt the wheelchair slightly backwards and move the wheelchair far enough onto the single carriageway so that the front wheels and footrests can be lowered without obstructing the kerb.

#### 4.3.5 Driving onto the kerb

Note: This operation is best performed with the assistance of a carer, and the best method is to drive to the kerb with the front of the carer, as this allows the direction of movement to be maintained and is less physically demanding for the carer.

##### 4.3.5.1 Stepping onto the kerb with the help of the attendant

- When the footrests are almost in contact with the kerb the pushchair should be inclined. The pram should be inclined enough so that the front wheels of the pram roll freely onto the pavement.
- Push the pushchair until the rear wheels make contact with the kerb. Then, while continuing to move the pushchair, lift the push handles while pushing the pushchair forward.

##### 4.3.5.2 Riding up onto the kerb in reverse with the help of the attendant

Note: For safety reasons, we do not recommend using this method.

#### 4.3.6 Ascending and descending stairs

Note: Avoid negotiating steps, particularly if no additional assistance is available. Use ramps or lifts if available.

Two carers are required for this activity, each capable of carrying at least half the combined weight of the wheelchair user and the wheelchair.

Each helper should be bent at the knees and holding the main frame of the wheelchair: the front helper should keep their hands under the front of the seat, while the rear helper should hold the frame of the wheelchair under the push handles.

#### **WARNING:**

Do not attempt to lift the wheelchair by its moving parts. Lifting by the moving parts of the trolley may cause injury to the user or damage to the trolley itself.

#### **THE LIFTING ACTION SHOULD BE CARRIED OUT WITH THE BACK STRAIGHT.**

Align the legs in relation to the wheelchair. Carefully lift it up or down. Be careful of the way your feet are placed - do not cross your feet.

#### 4.4 Transport

The Ergo Lite trolley (KM-2501 & KM-2512) has been accident tested in accordance with ISO7176-19 (Protocol No. MBK 09/0964) and is intended for use in vehicles that have been adapted for this purpose. The trolley must be securely fastened as described in its user manual. The Ergo Lite trolley (2501 & 2512) has been tested with a four-point belt attachment system. For further information on securing the pram, please refer to the manufacturer's instruction manual for the securing system. Karma always recommends that the user is transported on the seat of the vehicle secured with the seat belt.

#### **SUGGESTION:**

If possible, move the user to a seat in the vehicle.

##### 4.4.1 Before use

- Disconnect all detachable trolley accessories and store them in a safe place where they cannot move, so as not to create an additional hazard.

- The four-point attachment system of the trolley should adequately secure the lower frame of the trolley (hook attachment point ) to the transport vehicle.
- A restraint system using a harness must also be secured to the transport vehicle to ensure the safety of the user. Please refer to the system's instruction manual.
- The carer/person responsible for transport must check all fixing points. The carer/person responsible for transport must be fully familiar with the operation of the system and the weight limits for the system.
- See <http://www.unwin-safety.com/> for more information.

**WARNING:**

When transporting or lifting the wheelchair, do not hold it by moving parts or accessories (front/rear wheels, handles, front wheels, footrests, armrests and anti-tipping devices, etc.). The lifting procedure should be carried out on an upright back. Align the legs next to the wheelchair. Carefully lift the wheelchair up or down. Be careful of the way your feet are placed - do not cross your feet.

**4.4.2 Attachment method**

The Ergo Lite trolley (KM-2501 & KM-2512) has been tested with a four-point attachment system. Please refer to the system's instruction manual for further information on its operation.

**Front belt anchorage**

- Connect the straps around the frame of the wheelchair and the front wheel clamp as shown in the picture below. The attachment point is marked in yellow.
- Release the brakes and withdraw the trolley to tighten the front straps. Apply the trolley brakes (brakes should remain on during transport).

**Rear belt anchorage**

- Fit the back straps to the frame approximately 300 mm from the back of the wheelchair.
- Attach the hooks to the rear frame as shown in the diagram below. Tighten the straps until the trolley is securely restrained.


Front belt attachment

Rear belt attachment

**Seat belt attachment**

- Attach a 3-point safety belt for the user.
- The wheelchair user should be positioned facing forward and the wheelchair brake engaged.
- The pelvic belts should make full contact with the front of the body near the junction of the thigh and pelvis, as shown in the figure below (Correct seat belt position).
- The seat belt should not be held away from the body by elements of the wheelchair, such as armrests or wheels, as shown in the figure below (Incorrect seat belt position).

Correct seat belt position

Incorrect seat belt position

#### 4.5 safety belt

##### Seat belt installation and removal instructions

- Loosen the fourth and fifth bolts on each side of the seat upholstery and insert the seat belt between the fourth and fifth holes (Fig. 4.23).
- Keep the seat belt at an angle of 45° for comfort (Fig. 4.24).
- Slide the seat upholstery backwards and tighten the screws.


Fig. 23


Fig. 24

##### Seat belt length adjustment

- Determine the required lane length (Fig. 4.25).
- Slide the two parts of the belt through the buckle to lengthen or shorten the length of the belt (Fig. 4.26).
- Pull the two parts of the belt (Fig. 4.27).
- Slide the buckle into the plastic clamp so that the belt does not slip. Check that there is no space between the clamp and the buckle (Fig. 4.27).

##### SUGGESTION:

We recommend that the belt is fitted fairly tightly around the pelvis. The belt should be tightened in such a way that a finger can be inserted between the wearer's body and the belt. Check that the belt is fitted correctly, as recommended in the instructions. If this is not possible, please contact your supplier.

Fig. 4.25

Fig. 4.26

Fig. 4.27

Fig. 4.28

Belt maintenance:

- Once the user is placed in the wheelchair, the seat belt must be fastened.
- Adjust the belt to fit the wearer's body for comfort.
- If the belt is not used, it should be fastened so that it does not block the rear wheels when the pushchair is being driven.

Maintenance

The safety belt is made of 100% polyester and the buckle is made of plastic. The belt requires no maintenance, but should be inspected monthly for edge wear and weekly to check that it fits the wearer. If signs of wear are found, the belt should be replaced immediately.

## 5 INSPECTION AND MAINTENANCE

### 5.1 General inspection

Before driving, check the following items. If any malfunctions are detected, contact your KARMA dealer for an additional inspection of the wheelchair before use.

- Check brakes for wear, damage. Replacement or adjustment of the brakes should be carried out by an authorised technician.
- Check that the hand grips are not damaged and securely attached.
- Check that anti-tipping devices are not damaged and securely fitted.
- Before using the trolley, check that all bolts and nuts are tightened. Check all parts for wear and replace if necessary. Check the alignment of all parts.
- Check the brakes to make sure they are working properly.
- Wheels and tyres should be checked periodically for cracks and wear and replaced if necessary by an authorised technician.

#### **WARNING:**

Contact your Karma dealer to inspect the trolley if any faults are found.

### 5.2 Regular maintenance

In order to keep the pram in good condition, you should have it regularly inspected by an authorised karma dealer every six months from the date of purchase.

Maintenance/service charges may apply.

#### **SUGGESTION:**

If the trolley will not be used for a long period, it should be inspected annually.

Position	Every week	Every month
<b>GENERAL</b>		
Straight rollers (no excessive movement or pulling to one side)		V
<b>FRAME AND CROSSBARS</b>		V
Check the structure for small cracks/gaps Check bent frame or cross bars		V
<b>BRAKES</b>		
Check that they do not come into contact with the wheels when driving		V
Pivot points free from wear and play		V
switch on easily	V	

Prevent trolley movement	V	
<b>SEAT AND BACKREST</b>		
Check that they are not torn or sagging	V	
Check for play and damaged equipment		V
<b>FRONT WHEELS</b>		
Check that the front wheels turn smoothly		V
Check front wheels for worn bearings/mountings and possible tyre damage		V
<b>REAR WHEELS</b>		
That the quick-release axles close properly; lubricate if necessary.	V	
No sliding or pulling to one side		V
Check spokes for cracks, bends or fractures	V	
<b>ARMRESTS</b>		
Check armrest button for correct operation		V
Check that the armrests detach correctly and return to their position		V
<b>FOOTSTOOLS</b>		
Check that the footrest height adjustment screws are securely tightened.		V
Check that the footrest release lever locks correctly		V
<b>TYRES</b>		
Check for cracks, wear and tear	V	
For inflatable tyres, check that they inflate properly	V	
Check that the tyre tread depth is greater than 1 mm		V
Check tyre pressure		

**WARNING:**

As with any wheelchair, wheels and tyres must be inspected periodically for cracks and wear and replaced if necessary. Failure to maintain the correct tyre pressure significantly reduces brake performance.

**WARNING:**

Do not modify or disassemble the frame structure of the wheelchair. Doing so may cause injury. This also voids the warranty.

**6. OPTIONS AND ACCESSORIES**

Pillar I.V.  
Height adjustable from 135 cm to 200

Backrest bag  
Standard size  
Fits all Karma hand trucks

cm Fits all Karma hand trucks	
 <p>Drinks holder Adjusts to the size of the container (6 - 10 cm) Fits all Karma hand trucks</p>	 <p>Anti-tipping device Fits all Karma hand trucks</p>
 <p>Multifunctional seat belt Fits all Karma hand trucks</p>	

**WARNING:**

Fitting accessories to the wheelchair can disturb the overall stability and centre of gravity of the wheelchair. Thus, these operations must be carried out by an authorised technician.

**7. RECTIFICATION OF DEFECTS**

General trolley problems

You can check the following items yourself before returning the wheelchair for repair.

<b>Problem</b>	<b>Inspection sites and remedial action</b>
1. trolley changes direction right/left	<ul style="list-style-type: none"> <li>• Check tyres for correct pressure</li> <li>• Check front wheel clamping angle</li> <li>• Check that the rear wheels are mounted at an equal distance from the seat frame</li> </ul>
2. making slow turns	<ul style="list-style-type: none"> <li>• Check tyres for correct pressure</li> <li>• Check for loose nuts and bolts</li> </ul>
3. the front wheel spins	<ul style="list-style-type: none"> <li>• Check tyres for correct pressure</li> <li>• Check for loose nuts and bolts</li> <li>• Check front wheel clamping angle</li> </ul>
4. squeaks and crackles	<ul style="list-style-type: none"> <li>• Check for loose nuts and bolts</li> </ul>
5. trolley clearances	<ul style="list-style-type: none"> <li>• Check for loose nuts and bolts</li> </ul>

If you cannot solve the problem yourself, contact your Karma dealer

**8. SPECIFICATION**

Model	<b>Ergo-Lite(KM-2501)</b>
-------	---------------------------

Overall length (mm)	940 940	
Overall width (mm)	570 621	
Overall height (mm)	910	
Gross weight (kg) (pneumatic tyres)	8.3	8.4
Seat width (mm)	406(16")	457(18")
Seat depth (mm)	420	
Seat height (mm)	482	
Backrest height (mm)	420	
Height of armrest (mm)	215	
Calf length (mm)	350-400	
Length after folding (mm)	850	
Folded width (mm)	275	
Height when folded (mm)	710	
Maximum user weight (kg)	100 100	

Model	Ergo-Lite(KM-2512)			
	20"		14"	
Overall length (mm)	930	930	930	930
Overall width (mm)	605	656	570	621
Overall height (mm)	910	910	910	910
Gross weight (kg) (pneumatic tyres)	10.23	10,34	8,75	8,86
Seat width (mm)	406 (16")	457(18")	406 (16")	457(18")
Seat depth (mm)	420	420	420	420
Seat height (mm)	482	482	482	482
Backrest height (mm)	420	420	420	420
Height of armrest (mm)	215	215	215	215
Calf length (mm)	350	350	350	350
Length after folding (mm)	690	690	690	690
Width when folded (mm)	340	340	305	305
Height when folded (mm)	710	710	710	710
Maximum user weight (kg)	100	100	100	100

Available seat widths: 41 cm, 46 cm. Due to the product specifications, for the above wheelchair widths the actual seat width is respectively 37 cm and 43 cm.

Stability (according to ISO 7176-1:1999)			
100 kg dummy	Ergo Lite(KM-2501)	Ergo Lite(KM-2512)	
Forwards	16°	17.9°	18.1
Rearward - rear wheels locked	10.5°	10.1°	10.8°
Rearward - rear wheels unlocked	15°	15.5°	14.1°
Tipping prevention	11.5°	11.4°	11.9°
Side	16°	14.3°	13.3°

**WARNING:**

It is not permitted to drive the trolley across a hill (in either direction) with a gradient of more than 7 degrees.

## 9. CLEANING AND MAINTENANCE

### 9.1 CLEANING

To prolong the life of the trolley, it is recommended that it is cleaned periodically (preferably daily) especially if it is used in rainy or snowy weather. We recommend that you do not use the pushchair in rain or snow. If rain is accidentally encountered while moving the trolley, it is recommended to find a shelter. If it is not possible to find shelter, head to your destination.

#### 9.1.1 Frame

Wipe off with a damp cloth.

#### 9.1.2 Upholstery

- Use warm, soapy water to clean the upholstery.
- Upholstery cleaner/foam can also be used, but avoid solvents.

### 9.2 END-OF-LIFE RECYCLING

The end-of-life recycling plan should comply with environmental regulations. The items listed below are provided for information as local and national regulations may vary/change.

Material	Name of part	Method
Aluminium alloy	Frame	Recycling
Plastic	Modular elements	Recycling
Paper	Operating instructions	Recycling

- In many countries, there are a number of organisations that can be found online that accept used goods for charity. Please find such information and support those who have been a little less fortunate and protect the planet for future generations. If you do not have access to the Internet, you can use telephone information, ask your operator for help or get more information from your local information centre. Earth Day takes place every day.
- Contact your local recycling centre for information on battery recycling, as regulations may vary from country to country/region.

### NOTES:

In the event of a product-related 'serious incident' that directly or indirectly led to, could have led to or is likely to lead to any of the following:

- (a) the death of a patient, user or other person; or
- (b) temporary or permanent deterioration of the health of a patient, user or other person, or
- (c) a serious threat to public health

this 'serious incident' should be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is resident. For Poland, the competent authority is the Office for Registration of Medicinal Products, Medical Devices and Biocidal Products.

### NOTES:

If you experience pain, allergic reactions or other distressing, unclear symptoms related to the use of the medical device, consult a healthcare professional.

### HOW TO DISPOSE OF THE PRODUCT AT END-OF-LIFE

Once the device has been decommissioned, the medical device can be disposed of as normal municipal waste with the exception of electrical products - proceed in a manner appropriate to the disposal of electrical and electronic equipment.

### INFORMATION

This symbol indicates the maximum weight of the user

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### **Wytwórca:**

Karma Medical Products Co. Ltd. No. 2363 Sec. 2  
Da Shiue Road, Min-Hsiung Shiang, Chia-Yi 621, Tajwan

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Karma Medical Products (UK) Ltd., 28 Newfield Road, Hagley, Stourbridge, West Midlands, DY9 0JR, Wielka Brytania

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Model: .....  
Serial number: .....

Point of sale stamp

and a legible signature of the Seller: .....

1. ANTAR MEDICAL Sp. z o.o. 03-068 Warsaw, ul. Zawisłańska 43, hereby grants a 12-month warranty on the product, starting from the date of purchase by the Buyer.
2. During the warranty period ANTAR undertakes to carry out the repairs necessary to restore proper functioning of the product, within 14 days from the date of receipt of the product by the service centre, together with the proof of purchase and the warranty card.
3. ANTAR reserves the right to replace the product if the repair is deemed uneconomic. Repair or replacement of faulty products does not extend the warranty period.
4. If the product is found to be malfunctioning, please contact your retailer/distributor or (if it is impossible) send it to ANTAR (as specified in point 1).
5. Before use, read the instructions carefully and follow them. If the product is used improperly and contrary to recommendations, the warranty will be irrevocably lost.
6. The warranty does not cover damage caused by external mechanical force, load exceeding the maximum allowed weight of the user, as described in the manual, as well as contact of the product with alcohol, greasy substances or petrol. The warranty does not cover parts which are subject to normal wear or tear during use, such as the bearing, rubber parts and others.
7. The product delivered to the service in connection with the complaint cannot cause an epidemic hazard.
8. THIS WARRANTY IS ONLY VALID TOGETHER WITH THE PROOF OF PURCHASE (BILL, RECEIPT, INVOICE). IN CASE OF A COMPLAINT, A PROOF OF PURCHASE MUST BE ATTACHED.
9. Complaints submitted without the proof of purchase and a warranty card with the product serial number will not be considered.

**Manufacturer:** Karma Medical Products Co. Ltd.  
No. 2363 Sec. 2 Da Shiu Road,  
Min-Hsiung Shiang, Chia-Yi 621, Taiwan

**Authorised Representative:** Karma Medical Products (UK) Ltd, 28 Newfield Road, Hagley, Stourbridge, West Midlands, DY9 0JR, United Kingdom

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Výrobca: Karma Medical Products Co. Ltd.
Karma Medical Products (UK) Ltd, 28 Newfield Road, Hagley, Stourbridge, West Midlands, DY9 0JR, Spojené kráľovstvo